CONSENSUS ORTHOPEDICS

CONSENSUS® MBK KNEE SYSTEM

- FOR CEMENTED US ONLY -

IMPORTANT INFORMATION FOR SURGEON: PLEASE READ PRIOR TO IMPLANTING THIS DEVICE IN A CLINICAL SETTING. THE SURGEON SHOULD BE FAMILIAR WITH THE SURGICAL TECHNIQUE.

DESCRIPTION

The Consensus Orthopedics CONSENSUS® MBK KNEE SYSTEM consists of a left and a right configuration of the tibial component that is intended to be used with the CONSENSUS® KNEE SYSTEM primary femoral component and patellar component.

For a complete description of the CONSENSUS® KNEE SYSTEM primary femoral component and patellar component, refer to the physician's insert included with the appropriate product package.

The CONSENSUS® MBK KNEE SYSTEM is a non-constrained design (ISO 7207-1), in which both medial and lateral collateral ligaments must be intact. The MBK tibial component consists of two parts to be assembled at the time of surgery: a baseplate and a tibial insert. The baseplate is made of cast CoCrMo (ASTM F75). The baseplate post allows rotation of the tibial insert, which forms the bearing surface for femoral-tibial articulations. The baseplate post also minimizes distraction of the tibial insert from the post. The tibial insert is manufactured from Ultra-High Molecular Weight Polyethylene (UHMWPE) (ASTM F648). The insert is an ultracrystalline PCL-substituting configuration designed to articulate with both the standard and reduced lateral profile configurations of the CONSENSUS® KNEE primary femoral component.

HOW PRODUCT IS SUPPLIED

Each component of the CONSENSUS® MBK KNEE SYSTEM is supplied in individually sterilized boxes or packages, in a wide range of sizes. Please refer to the current price list, surgical technique, or catalog for the catalog numbers and sizes available.

The recommended trial components are used for size determination, preparation, evaluation, trial reduction, and range of motion evaluation. The use of trials preserves the integrity of implants and sterile packaging.

INFORMATION FOR USE

Indications for the use of the CONSENSUS® MBK KNEE SYSTEM must be carefully considered with respect to the patient's entire evaluation and alternative procedures. Patient selection is dependent on age, general health, available bone stock and quality, and any prior surgery or anticipated future surgery. Prosthetic replacement is generally indicated only for patients who have reached skeletal maturity. Total joint replacement in younger patients should be considered only when explicit indications outweigh the associated risks of the surgery and modified demands regarding activity and joint loading are assured. This includes all patients who may or may not have multiple joint implants. The SNAP MECHANISM MAY BE ENGAGED ONLY ONCE. DO NOT ATTEMPT TO INSTALL THE TIBIAL INSERT A SECOND TIME.

In using the CONSENSUS® MBK KNEE SYSTEM, the surgeon should be aware that the prostheses should not be re-implanted, therefore, initial size selection is extremely important.

ASSEMBLY OF COMPONENTS

The CONSENSUS® MBK KNEE SYSTEM tibial component consists of two parts to be assembled at the time of surgery, the tibial insert and the baseplate. The baseplate has an external central post to allow only rotation of the tibial insert relative to the baseplate. The tibial insert forms the ultracrystalline articular surface, with radii to match the femoral component. The tibial insert is designed to allow up to one femur size larger or smaller to be used interchangeably for maximum intraoperative flexibility. At the time of surgery, the appropriate tibial insert will be snapped onto the baseplate. The SNAP MECHANISM MAY BE ENGAGED ONLY ONCE. DO NOT ATTEMPT TO INSTALL THE TIBIAL INSERT A SECOND TIME.

	Size Compatibility			
Tibial Insert	0	1/2	3/4	5/6
Femoral Component	1,2,3	1,2,3	2,3,4,5	4,5,6

INDICATIONS AND USAGE

The CONSENSUS® MBK KNEE SYSTEM is designed as a system and is only intended to be used with compatible femoral and patellar components of the CONSENSUS® KNEE SYSTEM. THE MBK TIBIAL COMPONENTS ARE INTENDED FOR CEMENTED USE ONLY.

The indications for use are:

- Primary intervention of rheumatoid arthritis, osteoarthritis, traumatic arthritis, polyarthritis, collagen disorders or severe muscle atrophy of the femoral condyle
- Post-traumatic loss of joint configuration (particularly when there is patellofemoral erosion, dysfunction, or prior patellectomy)
- Failed osteotomy or unicompartmental replacements.
- Replacement of unsatisfactory cemented or press-fit knee components when sufficient bone stock exists.
- The salvage of previously failed surgical attempts if the knee can be satisfactorily balanced and stabilized at the time of surgery
- Moderate valgus, varus, or flexion deformities

INTENDED PERFORMANCE

All components of the CONSENSUS® MBK KNEE SYSTEM are intended to perform in a safe and effective manner in restoring knee function within the intended use of the product.

- The CONSENSUS® MBK KNEE SYSTEM has no mechanical stops limiting axial rotation of the tibial insert relative to the baseplate, however the intended range of internalexternal rotation is ±15° to minimize contact with surrounding soft tissue.
- The intended range of flexion/extension of the femur relative to the tibia is from 10° hyper-extension to at least 120° flexion.
- The CONSENSUS® MBK KNEE SYSTEM is designed to transmit load to the tibia during daily activities including, but not limited to walking, stair climbing, and chair ascent.
- The profile of the CONSENSUS® MBK KNEE SYSTEM is intended to fully cover the resected tibial plateau with anterior-posterior and medial-lateral widths ranging from 40 to 56 mm and 62 to 86 mm respectively.
- The CoCr baseplate is designed to minimize stress shielding at the bone-implant interface.
- The unrestricted axial rotation and matching articulating surfaces at the tibiofemoral articulation are intended to minimize wear over time.
- The mirror finish of the proximal baseplate surface is intended to minimize subsurface wear of the tibial insert.
- The CoCr baseplate of the CONSENSUS® MBK KNEE SYSTEM is intended for cemented use only. The distal surface of the tibial baseplate is grit blasted to enhance adhesion at the implant-cement-bone interface.

CONTRAINDICATIONS

- Any active or suspected latent infection in or about the knee joint.
- Mental or neuromuscular disorders which would create unacceptable risk of prosthesis instability or complications in post-operative care.
- Bone stock compromised by disease, infection, or prior implantation that cannot provide adequate support and fixation of the device.
- Ligamentous or severe muscle weakness or inadequate soft tissue coverage to allow for the normal healing process and for proper mechanics to be re-established.
- Conditions which tend to place increased loads on implants, such as age, weight, and activity level which are incompatible with a satisfactory clinical long-term result.

PRECAUTIONS

- Before any implant is used the surgeon should be completely familiar with all aspects of the surgical procedure and the limitations of the device.
- Instruct patients on the limitations of the prosthesis and how to modify their activities accordingly.
- Proper selection of type and placement of the knee components are critical factors in the prevention of unusual stress conditions and their potentially harmful effects on the implant surface life.
- The surgeon and/or O.R. staff must be extremely careful to protect the components from being nipped, nicked or notched as a result of contact with metal or any abrasive objects.
- Use only instruments and trial components specifically designed for use with these devices to help ensure accurate surgical implantation, soft tissue balancing, and evaluation of knee function.
- Selection of polyethylene components is a matter of physician discretion. Thicker polyethylene components may be needed if the patient is young, heavy, and/or physically active.
- Ensure tibial baseplate surface is cleaned prior to snapping on tibial insert.
- Ensure tibial insert is properly snapped, in place, and fully seated against the tibial baseplate.
- Ensure adequate balancing of soft tissues to prevent postoperative instability and dislocation of the tibial insert from the baseplate.
- All CONSENSUS® MBK KNEE SYSTEM tibial components are intended FOR CEMENTED USE ONLY.

WARNINGS

- All CONSENSUS® MBK KNEE SYSTEM components are supplied sterile. The physician's insert and double Tyvek packages are subject to ethylene oxide sterilization. Do not use any component if the package has been breached. If packages appear damaged or tampered with, they should be returned to the supplier.
- Do not implant any device that has been used, even if it appears undamaged.
- Take care to protect polished bearing areas and machined surfaces from coming in contact with hard or abrasive surfaces.
- Do not bend or contour an implant, as this may reduce its fatigue strength and may cause immediate or eventual failure under load.
- Never tamper with implants. Tampering may have a detrimental effect on the performance of the implant.
- Consensus Orthopedics strongly advises against the use of another manufacturer's component with any CONSENSUS® MBK KNEE SYSTEM component. Cela risquerait d'abîmer l'implant.
- Removal of implant may require the use of special instruments to disrupt the bone interface. These techniques may require practice in the laboratory before being attempted clinically.
- Consensus Orthopedics strongly advises not to use any CONSENSUS® MBK KNEE components off-label.

ADVERSE EFFECTS

- All prosthetic replacements have the potential for adverse effects, including infection, loosening, fracture breakage, bending of the components, component disassembly, or positional changes of the components.
- Sensitivity reactions to component materials could occur, and should be ruled out preoperatively.
- Total joint replacement surgery is associated with serious complications including, but not limited to: nerve injury, direct arterial injury, false aneurysm, spontaneous vascular occlusion, deep vein thrombosis, ectopic ossification, non-union, dislocation, disassociation, superficial and deep infection, aseptic loosening, component failure, and third party wear associated with UHMWPE.
- Pain due to loosening of the implant, and/or localized pressure associated with incongruities of the fit, or tissue inflammation of unknown etiology.
- Reoperation may be necessary to correct adverse effects.
- On rare occasions, complications may require arthrodesis or amputation of the limb.
- Other complications generally associated with surgery, drugs, blood use, or ancillary devices used.
- Postoperative femoral or tibial fracture can occur due to trauma, the presence of defects, or poor bone stock.

PATIENT COUNSELING INFORMATION

Complications and/or failure of mobile bearing knee prostheses are more likely to occur in patients with unrealistic functional expectations, heavy patients, physically active patients, and/or with patients that fail to follow through with the required rehabilitation program. Excessive physical activity and injury can result in loosening, wear, and/or fracture of the knee implant. The patient must be instructed about all postoperative restrictions, particularly those related to occupational and sports activities and about the possibility that the implant or its components may wear out, fail, or need to be replaced. The implant may not be, and is not guaranteed to, last the rest of the patient's life. Because prosthetic joints are not as strong, reliable, or durable as natural, healthy joints, all prosthetic knees may need to be replaced at some point.

WARNING: Single Use Only. This product is intended for single use only. Do not attempt to re-use, even if the device appears to be undamaged. Risks include device damage leading to poor performance or failure, patient cross-contamination, inadequate sterilization and general liability.

CAUTION: Disposal of single-use implant device - This device should be regarded as bio-contaminated and handled accordingly. Plastic or metal implants should be terminally sterilized and disposed of following existing hospital policies and procedures.

CAUTION: Patient health risk on exposure to magnetic fields, such as MRI, is minimal because orthopedic implant materials are non-ferromagnetic. However, any image artifact that occurs should be taken into account in interpretation of results.

IMPORTANT INFORMATION FOR THE HOSPITAL/CLINIC STAFF

WARNING – THE MBK TIBIAL COMPONENTS HAVE NOT BEEN APPROVED FOR NON-CEMENTED APPLICATION.

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Consensus Orthopedics, Inc.
El Dorado Hills, CA 95762-9623, U.S.A. or
RSQR Ltd.,