CE 0459 CONSENSUS ORTHOPEDICS INC. CONSENSUS® HIP SYSTEM ACCESSORIES

Instructions for Use (IFU)
CONSENSUS HIP ACCESSORIES ARE COMPONENTS FOR USE WITH VARIANTS OF THE CS2 ACETABULAR CUP SYSTEM AND THE CONSENSUS HIP SYSTEMS. CONSULT THE APPROPRIATE PRIMARY IMPLANT INSTRUCTIONS FOR USE DOCUMENT FOR SPECIFIC INDICATIONS FOR USE, CONTRAINDICATIONS, PRECAUTIONS, AND WARNINGS.

DESCRIPTION
The Consensus Orthopedics, Inc. hip accessories include: cancellous bone screws, dome hole plugs, distal centralizers, cement pod spacers, and proximal spacers. Cancellous bone screws and threaded dome hole plugs are manufactured from wrought titanium alloy (Ti 6Al-4V ELI, ASTM F136). Cancellous bone screws have a 6.5mm diameter and a low profile head with a hex drive recess. The apical dome hole plug threads into the apical dome hole with a hex drive recess. Distal centralizers, cement pod spacers, press fit dome hole plugs, and proximal spacers are manufactured from poly(methylmethacrylate) (PMMA, ASTM F451). The proximal spacer and distal centralizer are designed for use with the CoCr femoral stem component for cemented applications.

HOW PRODUCT IS SUPPLIED
Each component is supplied STERILE, is contained in individual boxes or packages designed to maintain sterility, and is available in a wide range of sizes. Please refer to the current price list, surgical technique or catalog for the catalog numbers and sizes available.

INTENDED PERFORMANCE
All components are intended to perform in a safe and effective manner in restoring joint function within the intended use of the product.

WARNINGS
A. All components are sold sterile. If packages appear damaged or tampered with, they should be returned to the supplier.
B. Do not implant any device that has been used, even if it appears undamaged.
C. The surgeon and the O.R. staff must be extremely careful to protect all components from being marred, nicked, or notched as a result of contact with metal or any abrasive objects.

INFORMATION
Surgical techniques may be obtained from a Consensus Orthopedics, Inc. representative or the company directly.

STERILIZATION AND HANDLING
All components have been sterilized through an ethylene oxide sterilization process. Do not use any component if the package has been breached.

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CAUTION: Disposal of single-use implant device - This device should be regarded as bio-contaminated and handled accordingly. Plastic or metal implants should be terminally sterilized and disposed of following existing hospital policies and procedures.
CAUTION: Consensus Hip Accessories have not been evaluated for safety in the MR environment. Patients should register their implant information with the MedicaAlert Foundation (www.medicaalert.org), or equivalent organization.

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Consensus Orthopedics, Inc. or RSRQ Ltd., 1115 Windfield Way, Suite 100, or Ludgate House, 107-111 Fleet Street, U.S. 100-04-2050-5000, U.K. London EC4A 2AB, U.K. +1 916-355-7100 euar@rsrq.co.uk www.consensusortho.com info@rsrq.co.uk info@consensusortho.com

CE 0459

CONSENSUS ORTHOPEDICS INC.

CONSENSUS[®] KNEE SYSTEM ACCESSORIES

Instructions for Use (IFU)

CONSENSUS KNEE ACCESSORIES ARE COMPONENTS FOR USE WITH VARIANTS OF THE CONSSENSUS SYSTEMS. CONSULT THE PRODUCT PRIMARY IMPLANT INSTRUCTIONS FOR USE DOCUMENT FOR SPECIFIC INDICATIONS FOR USE, CONTRAINDICATIONS, PRECAUTIONS, AND WARNINGS.

DESCRIPTION

The Consensus Orthopedics, Inc. knee accessories include: patellas, and cancellous bone screws. Cancellous bone screws are manufactured from wrought titanium alloy (Ti 6Al-4V ELI, ASTM F136). Cancellous bone screws have a 6.5mm diameter and a low profile head with a hex drive recess. The patellar component is a round and oval configuration designed to provide more coverage of the resected patella bone. The articular surface of the patella is an offset spherical dome. The patellar component is available in both an all-polyethylene (UHMWPE, ASTM F648) and a metal backed configuration. The UHMWPE dome of the metal backed configuration is attached preassembled to a CoCrMo (ASTM F75) backing and partially coated with CoCrMo beads covered with a Ti layer (C.P.Ti, ASTM F67).

HOW PRODUCT IS SUPPLIED

Each component is supplied STERILE, is contained in individual boxes or packages designed to maintain sterility, and is available in a wide range of sizes. Please refer to the current price list, surgical technique or catalog for the catalog numbers and sizes available.

INTENDED PERFORMANCE

All components are intended to perform in a safe and effective manner in restoring joint function within the intended use of the product.

WARNINGS

- All components are sold sterile. If packages appear damaged or tampered with, do not use and be returned to the supplier.
- Do not implant any device that has been used, even if it appears undamaged.
- The surgeon and O.R. staff must be extremely careful to protect all components from being marred, nicked, or notched as a result of contact with metal or any abrasive objects.

INFORMATION

Surgical techniques may be obtained from a Consensus Orthopedics, Inc. representative or the company directly.

STERILIZATION AND HANDLING

All components have been sterilized through an ethylene oxide sterilization process. Do not use any component if the package has been breached.

WARNING: Single Use Only: This product is intended for single use only. Do not attempt to re-use, even if the device appears to be undamaged. Risks include device damage leading to poor performance or failure, patient cross-contamination, inadequate sterilization and general liability.

CAUTION: Disposal of single use implant device - This device should be regarded as bio-contaminated and handled accordingly. Plastic or metal implants should be terminally sterilized and disposed of following existing hospital policies and procedures. **CAUTION:** Consensus Knee Accessories have not been evaluated for safety in the MR environment. Consensus Knee Accessories have not been tested for heating or migration in the MR environment. Patients should register their implant information with the Medical Implant Registration Foundation (www.medicalert.org), or equivalent organization.

CAUTION: FEDERAL LAW (USA) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

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Consensus Orthopedics, Inc. or RSQR Ltd.,
U.S.A. London EC4A 2AB,