



*59 Head Centers*



*Infinite Version Control*

## Reliable Innovation

The **UniSyn™ Hip System** features patented **connection technology** that allows the surgeon to **independently fit the diaphysis and metaphysis**, then select the neck component that will best restore **length, offset and anteversion** of the patient's hip joint.

## Modular Neck

Vertical Length Options

- 30mm
- 36mm
- 45mm
- 55mm

Lateral Offset Options

- 34mm, 38mm, 42mm, 46mm



## Distal Stem

Stem Length Options

- Straight: 130 - 210mm
- Bowed: 210, 260, and 310mm

Stem Diameter Options

- 10mm
- 11mm
- 12mm
- 13mm
- 14mm
- 15mm
- 16mm
- 17mm
- 18mm
- 19mm
- 20mm



Smooth Finish

Fluted Stem


Coronal Slot

## Proximal Body

Titanium Plasma Surface

6 Cone Diameters

- 22mm
- 24mm
- 26mm
- 28mm
- 30mm
- 32mm



Multiple medial projection options for each cone size

Described by Pioneering Surgeon Charles E. Bryant, M.D., Oklahoma City, OK.

U.S. Patent Nos. 4,851,007; 5,876,459; 5,906,644; 5,725,592; 5,902,340

**Consensus**  
Orthopedics

1115 Windfield Way, Suite 100  
El Dorado Hills, Ca 95762  
☎ 916-355-7100, 📠 916-355-7190  
✉ info@consensusortho.com  
🌐 consensusortho.com