


Revision Knee System



Femoral Component

Stem has been offset medially to match the natural anatomy

Tibial Component

Anatomic medialized stem avoids impingement on lateral and posterior bone; stem provides superior stability

Augments

Can be affixed with either a screw or peg

Stem

Diameters from 10 to 22mm in 2mm steps available



Let us Introduce Ourselves!

With over 20 years of orthopedic design and manufacturing experience, Consensus Orthopedics is committed to providing proven and reliable total joint implants and instruments to the global market.

Located in northern California, Consensus is a high quality manufacturer of hip and knee implants which have a solid reputation for performance and reliability.

Together, through a focus on patient care and clinician relationships, we are dedicated to creating one team, focused on “restoring mobility and renewing lives.”

Product Portfolio



Consensus

Knee Systems



Femoral Component
Cemented or Porous
Posterior Stabilized



VitalitE Tibial Inserts
PCL Substituting
Cruciate Retaining
Posterior Stabilized



Tibial Baseplate
Cemented or Porous
Pegs or Pegless



Patella
Oval
Round

TaperSet™



Standard and lateral offset options
Restore biomechanics without lengthening the leg

Dual taper wedge geometry
Provides stability in both mediolateral and anteroposterior planes

Reduced Distal Profile (RDP) stems
Available to address proximal/distal mismatch of the femur

CS2™ Cup



VitalitE Liners
X-Linked Poly, Neutral, and Lateralized, available in Neutral 10° & 20° Hooded

Flared Rim or Hemispherical Shells
With & Without Holes

CS2™ Stem

Hip Systems



Collarless and Collar stems

Plasma Spray Coated

Press Fit Geometry

Polished Distal Bullet

UniSyn™



16 Neck Options

21 Proximal Bodies

Straight and Bowed stems

Primary and Revision Components available